

INHOUSE SERVICES

Test Code	Test	Method	Specimen type	Additive	Container type	Amt. of sample	Temp	TAT	Special requirement
24HR	24 HR. URINE PROTEIN	Pyrogallol Red Method	Urine	Boric acid/Thymol	Sterile screw capped plastic container	NA	A/R	3 Hrs	24 Hr. Urine Collection
ABG	ABG	Direct Measurement	Whole blood(Arterial)	Heparin	Syringe	2ml	With cool pack if time is beyond 15 mins.	10 Mins after receiving	Should be reach to lab and processed immediately
AEC	ABSOLUTE EOSINOPHIL COUNT	Impedance +calculated Symex	Whole blood	K2 EDTA	Lavender top	2ml/3 ml	A/R	3 Hrs	Random
CBC	ABSOLUTE LYMPHOCYTE COUNT	Impedance +calculated Symex	Whole blood	K2 EDTA	Lavender top	2ml/3 ml	A/R	3 Hrs	Random

CBC	ABSOLUTE NEUTROPHIL COUNT	Impedance +calculated Symex	Whole blood	K2 EDTA	Lavender top	2ml/3 ml	A/R	3 Hrs	Random
APTT	ACTIVATED PARTIAL THROMBOPLASTIN TIME(APTT)	liquicelin-E with CaCl ₂ /Clot Detection	Citrated plasma	3.2% buffered sodium citrate	Blue Top	2.7 ml	A/R	2 Hrs	Random
ADA	ADA	End Point Method	Serum/Any Fluid/Other	Gel /Clot Activator	Yellow/Red cap	3.0 ml	A/R	2 Hrs	Random
AFB	AFB STAIN(ZN STAIN)	ZN Stain	Any Body Fluid/Sputum/an y Other	nil	Sterile screw capped plastic container	1-2 mL	A/R	1day	Random
SALB	ALBUMIN	Bromocresol Green	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
ACR	ALBUMIN/CREATININE RATIO(ACR)	Calculated Parameter	Urine	NA	Sterile screw capped plastic container	10 mL	A/R	2 Hrs	Random
AGR	ALBUMIN/GLOBULIN RATIO	Calculated Parameter	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
SAP	ALKALINE PHOSPHATASE	PNPP,AMP Buffer	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
AFP	ALPHA FETO PROTEIN (AFP)	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
ST	ALT/SGPT	UV with P5P	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	3 Hr	Random

SAMY	AMYLASE	Amylopectin	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	2 Hrs	Random
SAM	AMMONIA	Bromophenol Blue	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	1 Hrs	Random
ACCP	ANTI CCP	Latex Enhanced Turbidimetric Assay	Serum	Gel/Clot Activator	Yellow/Red cap	3.0ml	A/R	3 Hrs	Random
ASO	ANTISTREPTOLYSIN O(ASO)	Slide agglutination	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	2 Hrs	Random
ana	ANTI NUCLEAR ANTIBODY(ANA)	ELISA	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	2 Days	Random
CUSAF	ASCITIC FLUID C/S	Aerobic Culture	Ascitic Fluid	nil	Sterile screw capped plastic container	5 mL	A/R	48-72Hrs	Random
ASCITI C	ASCITIC FLUID R/E	Biuret, Hexokinase/ Light Microscopy	Ascitic Fluid	Sterile Container	Sterile screw capped plastic container	NA	A/R	3 Hrs	Random
OT	AST/SGOT	UV with P5P	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	3 Hrs	Random
HCG	BETA HCG	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow /Red cap	3.0ml	A/R	3 Hrs	Random
hc0h	BICARBONATE(ECO2)	Enzymatic	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	2 Hrs	Random
B.T	BLEEDING TIME(BT)	Ivy Method	Whole blood	NA	NA	NA	NA	On site	Random

CUSB	BLOOD CULTURE SENSITIVITY(AEROBIC)	Aerobic Culture	Blood	nil	Sterile screw capped plastic container	3-5 mL	A/R	3 Days	Random
FB	BLOOD GLUCOSE (FASTING)	GOD-POD/Hexokinase	Plasma	Fluoride	Grey Cap	3.0 mL	A/R	3 Hrs	Fasting
PP	BLOOD GLUCOSE (POST PRANDIAL)	GOD-POD/Hexokinase	Plasma	Fluoride	Grey Cap	3.0 mL	A/R	3 Hrs	Post Prandial
RB	BLOOD GLUCOSE (RANDOM)	GOD-POD/Hexokinase	Plasma	Fluoride	Grey Cap	3.0 mL	A/R	3 Hrs	Random
BG	BLOOD SUGAR GLUCOMETER	Electrochemical Technology	Whole Blood	NA	One Drop	NA	NA	NA	Random
BU	BLOOD UREA/BUN	Urease	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	3 Hrs	Random
BMT	BONE MARROW ASPIRATION	Staining & Slide Review	PSIS/Sternum	nil	NA	NA	A/R	2 Days	Random
F007	CA-125	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
SCAL	CALCIUM	ArsenazoIII	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	3 Hrs	Random
CHIK	CHIKUNGUNYA IgM	Card	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	2 Hrs	Random
CHL	CHLORIDE	ISE Direct/ISE Indirect	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	2 Hrs	Random
LP	CHOL/HDL RATIO	Calculated Parameter	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
C.T	CLOTTING TIME(C.T)	Lee -White Method	Whole blood	NA	NA	NA	NA	On site	Random

CORT	CORTISOL(MORNING, EVENING)	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
C3	COMPLEMENT3(C3)	Immunoturbidometric	Serum	Gel/Clot Activator	Yellow /Red cap	3.0ml	A/R	3Hrs	Random
C4	COMPLEMENT4(C4)	Immunoturbidometric	Serum	Gel/Clot Activator	Yellow/Red cap	3.0ml	A/R	3 Hrs	Random
CBC	COMPLETE BLOOD COUNT(CBC)	SLS/ sysmex, Impedance	Whole Blood	K2 EDTA	Lavender top	2ml/3 ml	A/R	3 Hrs	Random
CK	CREATINE KINASE(CK)	Rosalki	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	2 Hrs	Random
CKMB	CREATINE KINASE(CK-MB)	Enzymatic method	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	2 Hrs	Random
SC	CREATININE	Alkaline Picrate/Enzymatic	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	3 Hrs	Random
CRP	CRP(QUANTITATIVE)	Particle Enhanced Turbidimetric Immunoassay	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 m	A/R	2 Hrs	Random
CUSC	CSF C/S	Aerobic Culture	CSF	nil	Sterile screw capped plastic container	NA	A/R	48 - 72Hrs	Random
CSF	CSF R/E	Biuret, Hexokinase/ Light Microscopy	CSF	Sterile Container	Sterile screw capped plastic container	NA	A/R	2 Hrs	Random
550	DENGUE NS1 ANTIGEN(RAPID CARD)	Card	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	2 Hrs	Random

DC	DIFFERENTIAL COUNT(TLC& DLC)	Impedance / sysmex	Whole Blood	K2 EDTA	Lavender top	2ml/3 ml	A/R	2 Hrs	Random
SB	DIRECT BILIRUBIN	Spectrophotometric	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
DI	D-DIMER	Turbidometry	Serum	3.2% buffered sodium citrate	Blue Cap	3.0 ml	A/R	2 Hrs	Random
HDL	DIRECT HDL	Direct Measure/MgCl2	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	3 Hrs	Random
LDLC	DIRECT LDL	Calculated /Direct measurement	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
ESR	ESR	Manual/Westergren	Whole blood	3.8% buffered sodium citrate	Black top	1.6	A/R	2 Hrs	Random
FPER	FNAC	Staining & Slide Review	Any Site	nil	NA	NA	A/R	2 Days	Random
FSH	FSH	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 ml	A/R	3 Hrs	Random
FO	FOLATE	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 ml	A/R	3 Hrs	Random
H060	FERRITIN	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
R090	FREE T3	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
R091	FREE T4	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	3 Hrs	Random

FUCL	FUNGAL CULTURE	Conventional culture, Microscopy	Any Body Fluid/Sputum/an y Other	nil	Sterile screw capped plastic container	5ml	A/R	7days	Random
GPD	G6-PD	Kinetic Method	Whole Blood	K2 EDTA	Lavender top	3.0 mL	A/R	3 Hrs	Random
GGT	GAMMA GLUTAMYL TRASPEPTIDASE(GGT)	L-y-glutamyl-p-nitroanilide	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 MI	A/R	2 Hrs	Random
GENEX	GeneXpert MTB	CBNAAT	Sputum/Any Fluid	Nil	Sterile screw capped plastic container	NA	NA	1day	Random
	GENERAL HEALTH CHECK UP(CBC,RFT,URINE R/E,LIPID PROFILE,FBS)	CBC,RFT,URINE R/E,LIPID PROFILE,FBS	Serum, Whole Blood,Urine	K2 EDTA,Ge l	Lavender top/Yellow/r ed top, Sterile screw capped plastic container	3ml- 5ml	A/R	3 Hrs	Random
SG	GLOBULIN	Calculated Parameter	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 MI	A/R	3 Hrs	Random
GTT	GLUCOSE TOLERANCE TEST	GOD-POD/Hexokinase	Plasma	Fluoride	Grey Cap	3.0 MI	A/R	4 Hrs	Random
GS	GRAM STAIN	Gram Stain	Any Body Fluid/Sputum/an y Other	nil	Sterile screw capped plastic container	1-2 MI	A/R	1day	Random
HAV	HAV IgM	Card	Serum	Clot activators	Yellow top/Red top	3.0 MI	A/R	2 Hrs	Random

HB	HB & HCT	SLS/ sysmex	Whole Blood	K2 EDTA	Lavender top	2ml/3 ml	A/R	2 Hrs	Random
HB A1C	HBA1C	Immunoturbidometry	Whole Blood	K2 EDTA	Lavender top	3.0 ML	A/R	2 Hrs	NA
HCV	HCV(RAPID CARD)	Card	Serum	Clot activators	Yellow top/Red top	3.0 ML	A/R	2 Hrs	Random
HBEAG	HEPATITIS B e-ANTIGEN	Card	Serum	Clot activators	Yellow top/Red top	3.0 ML	A/R	2 Hrs	Random
HBSES	HEPATITIS B(ENHANCED CHEMILUMINESCENCE)	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 ML	A/R	3 Hrs	Random
HS	HEPATITIS B(RAPID CARD)	Card	Serum	Clot activators	Yellow top/Red top	3.0 ML	A/R	3 Hrs	Random
HCVES	HEPATITIS C(ENHANCED CHEMILUMINESCENCE)	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 ML	A/R	3 Hrs	Random
HEV	HEV IgM	Card	Serum	Clot activators	Yellow top/Red top	3.0 ML	A/R	2 Hrs	Random
HISS	HISTOPATHOLOGY SMALL	Microscopy	Tissue	Formalin	Unbreakable plastic containers	NA	A/R	10days	Random
HISM	HISTOPATHOLGY MEDIUM	Microscopy	Tissue	Formalin	Unbreakable plastic container	NA	A/R	10days	Random
HISL	HISTOPATHOLGY LARGE	Microscopy	Tissue	Formalin	Unbreakable plastic container	NA	A/R	10days	Random
HIVES	HIV(ENHANCED CHEMILUMINESCENCE)	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 ML	A/R	3 Hrs	Random

HCVRN A	HCV RNA Quantitative	CBNAAT	Plasma	K2 EDTA	Lavender top	2ml/3 ml	A/R	1 Day	Random
HBV	HBV DNA Quantitative	CBNAAT	Serum	Clot activators	Yellow top/Red top	3.0 ML	A/R	1 Day	Random
HIV	HIV(RAPID CARD)	Card	Serum	Clot activators	Yellow top/Red top	3.0 mL	A/R	2 Hrs	Random
INK	INDIA INK FOR PNEUMOCOCCI	India Ink	CSF	nil	Sterile screw capped plastic container	NA	NA	2 Hrs	Random
SB	INDIRECT BILIRUBIN	Direct Measure	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	3 Hrs	Random
TIBCS	IRON	Pyridyl azo Dye	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	2 Hrs	Random
FLU	INFLUENZA PCR	CBNAAT	Nasal and throat swab	VTM	VTM	NA	Cool Pack	1 day	Random
IgE	IMMUNOGLOBULIN IgE	Immunoturbidometric	Serum	Gel/Clot Activator	Yellow /Red cap	3.0ml	A/R	3 Hrs	Random
RFT	KIDNEY FUNCTION TESTS	Urea, Creatinine, Electrolytes	Serum	Gel/Red Top	Yellow cap/Red Cap	3.0 mL	A/R	3 Hrs	Random
FS	KOH (FUNGAL SMEAR)	KOH mount	Any Body Fluid/Sputum/any Other	nil	Sterile screw capped plastic container	1-2 mL	A/R	1 Day	Random
LP	LDL/HDL RATIO	Calculated Parameter	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random

LDH	LACTATE DEHYDROGENASE (LDH)	Pyruvate to Lactate	Serum, fluid	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
LH	LH	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow /Red cap	3.0ml	A/R	3 Hrs	Random
SL	LIPASE	Enzymatic with colipase	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	2 Hrs	Random
LP	LIPID PROFILE	Total Cholesterol,Triglycerides,HDL ,LDL ,VLDL, Total Cholesterol HDL ratio	Serum	Gel/Red Top	Yellow cap/Red Cap	3.0 mL	A/R	3 Hrs	Fasting
LFT	LIVER FUNCTION TEST	Total, Indirect Direct Bilirubin,ALT, AST,ALP	Serum	Gel/Red Top	Yellow cap/Red Cap	3.0 mL	A/R	3 Hrs	Random
SMAG	MAGNESIUM	Formazan Dye	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	3 Hrs	Random
MART	MALARIAL ANTIGEN(RAPID CARD)	Card	Whole Blood	K2 EDTA	Lavender top	3.0 mL	A/R	2 Hrs	Random
MALC	MALIGNANT CELL CYTOLOGY	Staining & Slide Review	Any Body Fluid/any Other	nil	Sterile screw capped plastic container	1-2 ml	A/R	1 Day	Random
MCH	MCH	Computed/Sysmex	Whole blood	K2 EDTA	Lavender top	2ml/3 ml	A/R	2 Hrs	Random
MCHC	MCHC	Computed /Sysmex	Whole blood	K2 EDTA	Lavender top	2ml/3 ml	A/R	2 Hrs	Random
MCV	MCV	Derived/Sysmex	Whole blood	K2 EDTA	Lavender top	2ml/3 ml	A/R	2 Hrs	Random

MICRO	MICROALBUMINURIA	Pyrogallol Red Method	Urine	NA	Sterile Container	3.0 mL	A/R	3 Hrs	Random
OSMOS	OSMOLALITY SERUM	Calculation	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	2 Hrs	Random
OSMO	OSMOLALITY URINE	Calculation	Urine	nil	Sterile screw capped plastic container	10 mL	A/R	3 Hrs	Random
PAPS	PAP SMEAR EXAMINATION	Staining & Slide Review	Cervical smear(transforma tion zone)	nil	NA	NA	A/R	2 Days	Random
PCV	PCV	Computed Pulse Height Sysmex	Whole blood	K2 EDTA	Lavender top	2ml/3 ml	A/R	3 Hrs	Random
PBF	PERIPHERAL SMEAR	Blood film/light microscopy	Whole blood	K2 EDTA	Lavender top	2ml/3 ml	A/R	1 Day	Random
MP	PERIPHERAL SMEAR FOR MALARIAL PARASITE	Thick and thin film/ light microscopy	Whole blood	K2 EDTA	Lavender top	2ml/3 ml	A/R	3 Hrs	Random
SPH	PHOSPHORUS	Phosphomolybdate reduction	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	2 Hrs	Random
PLAT	PLATELET COUNT	Impedance / sysmex	Whole Blood	K2 EDTA	Lavender top	2ml/3 ml	A/R	3 Hrs	Random
CUSPF	PLEURAL FLUID C/S	Aerobic Culture	Pleural Fluid	nil	Sterile screw capped plastic container	5 mL	A/R	48-72 Hrs	Random
PF	PLEURAL FLUID R/E	Biuret, Hexokinase/ Light Microscopy	Pleural Fluid	Sterile Container	Sterile screw capped	NA	A/R	3 Hrs	Random

					plastic container				
KA	POTASSIUM	ISE Direct/ISE Indirect	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
PRO	PROLACTIN	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
NTBNP	PRO BNP	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 ml	A/R	3 Hrs	Random
PCT	PROCALCITONIN(PCT)	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
SPSA	PROSTATE SPECIFIC ANTIGEN TOTAL(PSA)	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
PCR	PROTEIN/CREATININE RATIO(PCR)	Calculated Parameter	Urine	NA	Sterile screw capped plastic container	10 mL	A/R	3 Hrs	Random
PTI	PROTHROMBIN TIME(PT/INR)	Thromborel S/Clot detection	Citrated plasma	3.2% buffered sodium citrate	Blue Top	2.7 ml	A/R	2 Hrs	Random
CUSP	PUS/TISSUE/OTHER C/S	Aerobic Culture	Pus/Tissue/Any Other	nil	Sterile screw capped plastic container	NA	A/R	48-72 Hrs	Random
RAT	RAPID TEST FOR COVID 19	RAPID	Nasal and throat swab	VTM	VTM	NA	Cool pack	30 mts	Random

RBC	RBC	Impedance sysmex	Whole blood	K2 EDTA	Lavender top	2ml/3 ml	A/R	3 Hrs	Random
RC	RETICULOCYTE COUNT	Supravital staining/Light microscopy	Whole blood	K2 EDTA	Lavender top	2ml/3 ml	A/R	1 Day	Random
RA	RHEUMATOID FACTOR(RA FACTOR)	Slide agglutination/Immunoturbidometry	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	3 Hrs	Random
ELCT	S. ELECTROLYTES	Sodium, Potassium,Chloride	Serum	Gel/Red Top	Yellow cap/Red Cap	3.0 mL	A/R	3 Hrs	Random
SARS	SARS CoV 2 RNA	CBNAAT	Throat and nasal swab	VTM medium	Special Vial	NA	A/R	1 Day	Random
SCRUB T	SCRUB TYPHUS	RAPID	Serum	Gel/Red Top	Yellow cap/Red Cap	3.0ml	A/R	2 Hrs	Random
SA	SEMEN ANALYSIS	Conventional method	Semen	nil	Dry, sterile screw capped, plastic container ,Condoms are not accepted.	NA	A/R	2 Hrs	3 -5days Abstinence
NA	SODIUM	ISE Direct/ ISE Indirect	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
CS	STOOL C/S	Aerobic Culture	Stool	nil	Sterile screw capped plastic container	1 gm	A/R	48-72 Hrs	Random
HDFC	STOOL FOR HANGING DROP	Light Microscopy	Stool	nil	Sterile screw capped	1 gm	A/R	1 Hr	Random

					plastic container				
OB	STOOL OCCULT BLOOD	Benzidine Test	Stool	nil	Screw capped plastic container	1gm	A/R	2 Hrs	Random
SRS	STOOL REDUCING SUBSTANCE	Benedict's	Stool	nil	Sterile screw capped plastic container	1gm	A/R	3 Hrs	Random
SRE	STOOL ROUTINE AND MICROSCOPIC	Visual and light microscopy	Stool	nil	Sterile screw capped plastic container	1gm	A/R	3 Hrs	Random
SYNOV	SYNOVIAL FLUID R/E	Biuret, Hexokinase/Light Microscopy	Synovial Fluid	Sterile Container	Sterile screw capped plastic container	NA	A/R	3Hrs	Random
Z480	TETRAIODOTHYRONINE(T4)	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	3 Hrs	Random
TEST	TESTOSTERONE TOTAL	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	3 Hrs	Random
R052	THYROID STIMULATING HORMONE(TSH)	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	3 Hrs	Random
F1	THYROID FREE(FT3,FT4,TSH)	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
Z480	THYROID TOTAL(T3,T4,TSH)	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random

TIBCS	TIBC	Chromazurol B	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	3 Hrs	Random
S196	TISSUE TRANSGLUTAMINASE IGA(TTG)	CLIA	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	2 days	Random
SB	TOTAL BILIRUBIN	Diphylline, Diazonium Salt	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	3 Hrs	Random
LDL	TOTAL CHOLESTEROL	Cholesterol oxidase,esterase,peroxidase	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
SPT	TOTAL PROTEIN	Biuret end point	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	3 Hrs	Random
25H	TOTAL VITAMIN D(25-OH)	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
TC	TOTAL WBC COUNT	Impedance / sysmex	Whole Blood	K2 EDTA	Lavender top	2ml/3 ml	A/R	3 Hrs	Random
STR	TRIGLYCERIDE	Enzymatic, End point	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	3 Hrs	Random
Z480	TRIIODOTHYRONINE(T3)	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random
TROPIE S	TROPONIN -I (ES)	Enhanced Chemiluminescence	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	30 mins	Random
TYP	TYPHI DOT	Card	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	2 Hrs	Random
SUA	URIC ACID	Uricase	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	3 Hrs	Random

PT	URINE PREGNANCY (HCG)	Reagent strip	Urine	nil	Sterile screw capped plastic container	10 mL	A/R	2 Hrs	First Morning sample preferred
UPH	URINE (PH)	Reagent strip double ph indicator method	Urine	nil	Sterile screw capped plastic container	10 mL	A/R	2Hrs	First Morning/Random sample
URE10	URINE 10 PARAMETERS(glucose,bilirubin,ketone body,specific gravity,pH,protein,urobilinogen,nitrate,leucocytes)	Reagent strip /Light microscopy	Urine	nil	Sterile screw capped plastic container	10 mL	A/R	2 Hrs	First Morning sample preferred
UBP	URINE BILE PIGMENTS	Reagent strip /Fouchet's Test	Urine	nil	Sterile screw capped plastic container	10 mL	A/R	3 Hrs	First Morning/Random sample
UBS	URINE BILE SALTS	Hay's Test	Urine	nil	Sterile screw capped plastic container	10 mL	A/R	3 Hrs	First Morning/Random sample
UCT	URINE CREATININE 24 Hrs.	Alkaline Picrate	Urine	Boric acid/Thymol	Sterile screw capped plastic container	NA	A/R	3 Hrs	24 Hr. Urine Collection
CRU	URINE CREATININE SPOT	Alkaline Picrate	Urine	nil	Sterile screw capped plastic container	NA	A/R	3 Hrs	Random

US	URINE GLUCOSE	Reagent strip/Benedict's	Urine	nil	Sterile screw capped plastic container	10 mL	A/R	2 Hrs	First Morning/Random sample
UACE	URINE KETONES	Reagent strip/ Rothera's	Urine	nil	Sterile screw capped plastic container	10 mL	A/R	2 Hrs	First Morning/Random sample
UM	URINE MICROSCOPIC EXAMINATION	Light microscopy	Urine	nil	Sterile screw capped plastic container	10 mL	A/R	2 Hrs	First Morning/Random sample
POR	URINE PORPHOBILINOGEN	Ehrlich's Reagent	Urine	nil	Sterile screw capped plastic container	10 mL	A/R	3 Hrs	First Morning/Random sample
UKA	URINE POTASSIUM (K)SPOT	Direct /Indirect ISE	Urine	nil	Sterile screw capped plastic container	NA	A/R	2 Hrs	Random
UK24	URINE POTASSIUM 24 Hrs.	Direct /Indirect ISE	Urine	Boric acid/Thymol	Sterile screw capped plastic container	NA	A/R	2 Hrs	24 Hr. Urine Collection
URE	URINE R/E	Reagent strip /Light microscopy	Urine	nil	Sterile screw capped plastic container	10 mL	A/R	2 Hrs	First Morning sample preferred

US	URINE REDUCING SUGAR	Benedict's	Urine	nil	Sterile screw capped plastic container	10 mL	A/R	2 Hrs	First Morning/Random sample
UNA	URINE SODIUM (Na)SPOT	Direct /Indirect ISE	Urine	nil	Sterile screw capped plastic container	NA	A/R	3 Hrs	Random
UN24	URINE SODIUM 24 Hrs.	Direct /Indirect ISE	Urine	Boric acid/Thymol	Sterile screw capped plastic container	NA	A/R	3 Hrs	24 Hr. Urine Collection
UU	URINE UROBILINOGEN	Reagent strip	Urine	nil	Sterile screw capped plastic container	10 mL	A/R	2 Hrs	First Morning/Random sample
UALB	URINE: PROTEIN	Reagent strip/SSA	Urine	nil	Sterile screw capped plastic container	10 mL	A/R	2 Hrs	First Morning/Random sample
USGR	URINE: SPECIFIC GRAVITY	Reagent strip	Urine	nil	Sterile screw capped plastic container	10 mL	A/R	2 Hrs	First Morning/Random sample
VDRL	VDRL	Card	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	2Hrs	Random
VB12	VITAMIN B12	Enhanced Chemiluminescence	Serum	Clot activators	Yellow top/Red top	3.0 mL	A/R	3 Hrs	Random

VM	VIRAL MARKER	Enhanced Chemiluminescence	Serum	Clot activators	Yellow top/Red top	3.0 mL	A/R	3Hrs	Random
LP	VLDL	Calculated Parameter	Serum	Gel/Clot Activator	Yellow/Red cap	3.0 mL	A/R	2 Hrs	Random
WI	WIDAL	Slide/tube agglutination	Serum	Gel/Clot Activator	Yellow /Red cap	3.0 mL	A/R	2 Hrs	Random